



Committed to customer
Committed to future...

Veterinary Mobile 100mA



Touch Screen
Technology with 16
Bit DSP Controller

100 mA Mobile X-Ray

Manufactured & Marketed by:
MediRay
Healthcare

Office Plant:
Link Road, Near Car Market,
Sector-28, Faridabad, Haryana 121008
Tel: +91 9999665300

Overseas Enquiry:
Mobile: +91 9990778887
Whats app: +91 9990778887
E-mail: vinod.chaudhary@mediray.in

As ISO 15001:2015 and ISO 13000 Certified Company

Mediray Healthcare Nepal Pvt. Ltd.
Katu No.251, Rapti Gaunpalika,
Ward No.2, Saptari, Province 2, Nepal
Tel: +977 9807740444

Mediray Europe:
Ul. Mostowa 1, 31-061
Krakow, POLAND
Tel: +48 71 076269

Web: www.mediray.in
E-mail: info@mediray.in






UltimaX SERIES Generator Specifications

NEW

NEW

Specifications	UltimaX- 100	UltimaX- RX100 (Single Tank)	UltimaX- RX-200 (Single Tank)
Power Rating	8 KW	10 KW	16 kW
Kvp Range	30 to 100 kvp (1kvp step)	30 to 100 kvp (1kvp step)	30 to 125 kvp (1kvp step)
Timer Range (Second)	0.01 to 10 sec (1000 steps)	0.01 to 10 sec (1000 steps)	0.01 to 10 sec (1000 steps)
Technical Selection	25mA, 50mA, 100mA	25mA, 50mA, 100mA	50mA, 100mA, 150mA, 200mA
Tube	BEL/Imported Stationary anode X-Ray tube	BEL/Imported Stationary anode X-Ray tube	Toshiba/BEL Rotating anode X-Ray tube
APR Program (Pre-Programmed & Customs)	648	648	648
Input line (50/60 Hz)	Single Phase (170-260 VAC)	Single Phase (170-260 VAC)	Single Phase (170-260 VAC)
Remote	Optional	Optional	Optional
Stored Energy (SE)	Optional	Optional	Optional

Due to continuous product improvements, specifications & appearance are subject to change without prior notice.



Manufactured & Marketed by :

MediRay
Healthcare

An ISO 13485:2003 and 9001:2008 Certified Company

Office/Plant :
Link Road , Near Car Market,
Sector-28, Faridabad, Haryana 121008
Tel.: +91 9599565300

Web: www.mediray.in
E-mail : info@mediray.in

