



Committed to customer
Committed to future...



Veterinary X-ray 10mA



Wall mount



Portable



Manufactured & Marketed by :

MediRay
Healthcare

An ISO 13485:2003 and 9001:2008 Certified Company

Office/Plant :

Link Road , Near Car Market,
Sector-28, Faridabad.Haryana 121008
Tel.: +91 9599565300

Mediray Healthcare Nepal Pvt. Ltd.
Kita No.251, Rupani Gaunpalika,
Ward No.2, Saptari, Province 2, Nepal
Tel.: +977 9807740444

Overseas Enquiry:

Mobile : +91-9990778887
Whats app: +91-9990778887
E-mail : vishnu.chaudhary@mediray.in

Mediray Europe:
Ul. Mostowa 1, 31-061
Krakow, POLAND
Tel.: +48 791076269

Web: www.mediray.in
E-mail : info@mediray.in



SALIENT FEATURES

- ⊕ Hand Held Light Weight Portable
- ⊕ Perfect positioning in your Hand
- ⊕ Soft touch timing selection.
- ⊕ Consecutive exposures
- ⊕ last exposure selection retained in memory.
- ⊕ Lead coated Radiography Long Cone.
- ⊕ RVG Mode for Digital Imaging
- ⊕ Digital Display.



TECHNICAL SPECIFICATIONS

X-Ray output (Digital)	: 55 K.V. 4mA
Tube filtration	: 1.5mm All
Focal Point	: 0.8mm
Focus skin distance	: 10cm
Safety leakage radiation	: 30mR/h at 1 m. distance
Time Range (Digital)	: 0.0 Sec. to 2.0 Sec.
Power Supply	: 230v. 50Hz. 5 Amp.

SALIENT FEATURES

- ⊕ Microprocessor Controlled Dental X-Ray
- ⊕ The X-Ray emitter is connected to Scissor Arm and is movable both in Vertical and Horizontal directions.
- ⊕ The rotation of the Tube Head about its axis in 360° in Horizontal and 270° in Vertical plane.
- ⊕ Anatomical Programming Radiography (30 programs) with Selection.
- ⊕ Option to select three patient size



TECHNICAL SPECIFICATIONS

X-Ray output (Digital)	: 60 K.V. 8mA
Tube filtration	: 1.5mm All
Focal Point	: 0.8mm
Focus skin distance	: 15cm
Safety leakage radiation	: 30mR/h at 1 m. dis.
Time Range (Digital)	: 0.0 Sec. to 3.0Set.
Line voltage compensation:	: 190v to 260v