

## VETERINARY X-RAY PORTABLE

- > Truly Portable Carry away box type Machine
- > Voltage Over Load Protection
- > Tape Built FFD Measurement
- > Micro Processor Based Digital with Accuracy of 0.013sec.
- > Exposure Time can be set by Knob
- > Hand Switch for exposure in Single Touch
- > Weight 20KG



**UltimaX-50**



**UltimaX-30**

Manufactured & Marketed by :

**MediRay**  
Healthcare

An ISO 13485:2003 and 9001:2008 Certified Company

**Office/Plant :**

Link Road , Near Car Market,  
Sector-28, Faridabad.Haryana 121008  
Tel.: +91 9599565300

**Mediray Healthcare Nepal Pvt. Ltd.**

Kita No.251, Rupani Gaunpalika,  
Ward No.2, Saptari, Province 2, Nepal  
Tel.: +977 9807740444

**Overseas Enquiry:**

Mobile : +91-9990778887  
Whats app: +91-9990778887  
E-mail : vishnu.chaudhary@mediray.in

**Mediray Europe:**

Ul. Mostowa 1, 31-061  
Krakow, POLAND  
Tel.: +48 791076269

Web: [www.mediray.in](http://www.mediray.in)

E-mail : [info@mediray.in](mailto:info@mediray.in)

All are ISO 13485:2003 Certified Company



Specifications	UltimaX 50	UltimaX 30
<b>Output</b>	50mA at 50 Kv	30mA at 55 Kv
	40mA at 65 Kv	20mA at 75 Kv
	35mA at 65 Kv	15mA at 85 Kv
	30mA at 85 Kv	
<b>Time</b>	Microprocessor	Microprocessor
	0.1 to 6 sec.	0.2 to 5 sec.
	60 Steps	60 Steps
<b>Tube</b>	1.44mm Focal Spot	1.44mm Focal Spot
	Imported	Imported
<b>Input Voltage</b>	230V , 5AMP	230V , 5AMP
<b>Compensation</b>	210 to 250V	210 to 250V
<b>BeamLimitation</b>	Light Beam Collimator	Light Beam Collimator
<b>Weight</b>	22KG	20KG

**Special Features :**

- Truly Portable Carry away box type Machine
- Voltage Over Load Protection
- Tape Built FFD Measurement
- Micro Processor Based Digital with Accuracy of 0.015Sec.
- Exposure Time can be set by Knob
- Full Proof Radiation Protection
- Line Voltage Compensation
- Hand Operated switch if provided for X-Ray exposure in Single Touch

**STANDARD ACCESSORIES :**

- Power Supply Cable, Exposure Switch with Cable, Beam Limitation Device, Carrying Case (Optional)

**STANDS**

- ALLO Ratche Stand : - Light Weight (Optional)
- Dismantable Stand Convenient
- Ratch Type Stand

