



**EW-100**  
Washer

# EW-100 Washer



## Technical Specifications

Plate type	96-well plate
Washing head	8/12 type
Residual aspiration volume	< 2ul
Dispensing volume	10~3000ul, step by 1ul
No. of washing cycle	1~99 adjustable
Soak time	0~999s
Shaking time	0~999s
Temperature	5°C~35°C
Humidity	15%~85%
Power supply	AC 110V, 60Hz AC 220V, 50Hz
Dimension (mm)	400(L) x 310(W) x 169(H)
Weight	8 kg



### Mediray Healthcare Pvt. Ltd.

An ISO 13485:2003 Certified Company

Regd. Off: 513, Pkt-2, OM Apartment, Sector 14, Dwarka, New Delhi-110075 (INDIA)

Corp Off: Plot No.20, Pocket-4, B-Block, Sector-23B, Near Sec-21 Metro Station, Dwarka, New Delhi-110075 (INDIA)

Factory: E-179, Trans Delhi Signature City, Sec-A-7, UPSIDC, Loni, Ghaziabad(UP) 201102 (INDIA)

web: [www.mediray.in](http://www.mediray.in)

### DISTRIBUTOR

